

FabGennix International, Inc.

9191 Kyser Way Bldg. 4 Suite 402 Frisco, TX 75033

Tel: (214)-387-8105, 1-800-786-1236 Fax: (214)-387-8105 Email: info@fabgennix.com Web: www.FabGennix.com

PDE2A1 Western Blot Positive Control

Catalog Number: PC-PD2A1

Lot Number:

General Information

Product	PDE2A1 Positive Control
Description	Western Blot Positive Control Phosphodiesterase 2A1
Accession #	Uniprot: O00408 NCBI: NP_001230713.1
Verified Applications	WB
Immunogen	Purified protein.
Alternative Nomenclature	Phosphodiesterase 2A2

Physical Properties

Quantity	5 Applications
Volume	Inquire
Form	Positive Control in SDS-sample buffer
Storage	Store at -20°C for long term storage.
Form	

Recommended Dilutions

Western Blot 1:250-1:1,000

Related Products Catalog

PD2A1 Antibody	PD2A1-101AP
BIOTIN-Conjugated	PD2A1-BIOTIN
FITC-Conjugated	PD2A1-FITC
Antigenic Blocking Peptide	P-PD2A1
Pan PDE2A Antibody	PD2A-101AP
PDE2A2 Antibody	PD2A2-201AP
PDE2A3 Antibody	PD2A3-301AP
PDE2A4 Antibody	PD2A4-121AP

Tel: (214)-387-8105, 1-800-786-1236 Fax: (214)-387-8105 Email: info@fabgennix.com Web: www.FabGennix.com

Overview:

Cyclic nucleotides are important intracellular second messengers which play important role variety of signal transduction process. The cyclic nucleotides are hydrolyzed and compartmentalized by a family of enzymes called phosphodieterases. One of the phosphodiesterases that compartmentalize and hydrolyze cAMP and cGMP into AMP and GMP respectively, are phosphodiesterase type 2. There are two members of the PDE2A gene that are cloned and designated as PDE2A and PDE2B (1). Each of the PDE2A and PDE2B genes are 2793 Bases and 930 amino acids with an apparent MW of 100 kDA. The catalytic domain of PDE2A is homologous to all other known PDE family members. The PDE2A has specificity for cAMP with Km of 2.4 µM and its activity is modulated by the presence of cGMP. In contrast, the PDE2B does not hydrolyze cGMP nor is its activity altered by the presence of cGMP. Western blot analysis of PDE2A reveals a variety of tissues including neocortex, cerebellum, heart, kidney lungs, pulmonary artery and skeletal muscle (2). The PDE2A expression was evident in venous arterial endothelial cells but not in arterial endothelial cells (2). The PDE2A expression was also noted in corpus cavernosum along with a wide repertoire of other PDEs enzymes (3). There are currently 4 isoforms of the PDE2A protein denoted as PDE2A1, PDE2A2, PDE2A3 and PDE2A4.

FabGennix provides antibodies to many PDE family members including PDE-selective, family subtype-selective and family-subtype-variant selective antibodies for the detailed analysis of cyclic nucleotide signaling pathways. PDE2A1-selective antibody has been generated against N-terminus amino acid sequence taken within the first 50 amino acids of human PDE2A1 protein. The PDE2A1-specific antiserum has no cross reactivity against PDE2B or with other PDE family members tested so far. Antigenic blocking peptides (P-PD2A1) and western blot positive controls (PC-PD2A1) are available. PD2A1 western blot positive control is prepared in ready-to-use SDS-sample buffer. Antibodies can be conjugated to secondary enzymes or fluorophores upon request at nominal costs. For a complete listing of all FabGennix products and services please visit http://fabgennix.com.

References

- 1. Rascon A., Soderling SH, Schaefer and Beavo JA. Proc. Natl. Acad. Sci. USA 2002; 299, 4714-4719.
- 2. Sadhu K., Hanseley, K, V. L. Florio, and Folda, SA. J. Histochem and Cytochem. 47, 895-906, 1999.
- 3. Kuthe A., Widenroth, A., Magret J. H., et. al., J. Urol. 2001. 165, 280-283.

* For users who may require large amounts of the products listed above, please inquire about bulk material discounts. This Product is for Research Use Only and is NOT intended for use in humans or clinical diagnosis.

Tel: (214)-387-8105, 1-800-786-1236 Fax: (214)-387-8105 Email: info@fabgennix.com Web: www.FabGennix.com