

FabGennix International, Inc. 9191 Kyser Way Bldg. 4 Suite 402 Frisco, TX 75033

Tel: (214)-387-8105, 1-800-786-1236 Fax: (214)-387-8105

Email: info@fabgennix.com Web: www.FabGennix.com

Rabbit Polyclonal Anti-PDE2A antibody

Catalog Number: PD2A-101AP

Lot Number:

General Information

PDE2A Antibody
Phosphodiesterase 2A, cGMP-stimulated
Antibody Affinity Purified C-epitope
Uniprot: UPI000003B340
GenBank: EAW74861.1
CM, ELISA, ICC, IF, IHC, IP, WB
Human, Mouse, Rat
Rabbit
Synthetic peptide corresponding to unique C-
terminal amino acid sequence on PDE2A gene.
This antibody will label PDE2A1, PDE2A2,
PDE2A3 and PDE2A4 variants.
5"-cyclic phosphodiesterase antibody, cGMP
dependent 3'5' cyclic phosphodiesterase antibody,
Phosphodiesterase 2A cGMP stimulated antibody

Physical Properties

Quantity	100 µg
Volume	200 µl
Form	Affinity Purified Immunoglobulins
Immunoglobulin & Concentration	0.7 mg/ml IgG in antibody stabilization buffer
Determinant	C-epitope
Storage	Store at -20°C for long term storage.

Recommended Dilutions


DOT Blot	1:10,000
ELISA	1:10,000
Immunocytochemistry	1:150-1:1,000
Immunofluorescence	1:150-1:1,000
Immunohistochemistry	1:250-1:1,000
Immunoprecipitation	1:200
Western Blot	1:250-1:1,000

Related Products

Catalog

FITC-Conjugated	PD2A.c-FITC
BIOTIN-Conjugated	PD2A.c-BIOTIN
PDE2A N-epitope	PD2A-112AP
Antigenic Blocking Peptide	P-PD2A.c
Western Blot Positive Control	PC-PD2A

Application Verification:



Immunostaining of PDE2A expressing cells with PD2A-101AP. 1:250 antibody dilution in DiluObuffer.

IF with PD2A-101AP with whole-mount mouse nose. Green: Grueneberg ganglion neurons. Red: antibody immunoreactivity. Staining performed with blocking peptide P-PD2A.n. Data provided by Dr. Kroos Lab Caltech, Pasadena, CA.

Dilutions are for reference only. Applications not listed above are not necessarily precluded from working with this antibody. Investigators intending to use an application that has not been verified can request a complimentary sample.

Overview:

Cyclic nucleotides are important intracellular second messengers which play important role variety of signal transduction process. The cyclic nucleotides are hydrolyzed and compartmentalized by a family of enzymes called phosphodieterases. One of the phosphodiesterases that compartmentalize and hydrolyze cAMP and cGMP into AMP and GMP respectively, are phosphodiesterase type 2 There are two members of the PDE2A gene that are cloned and designated as PDE2A and PDE2B (1). Each of the PDE2A and PDE2B genes are 2793 Bases and 930 amino acids with an apparent MW of 100 kDA. The catalytic domain of PDE2A is homologous to all other known PDE family members. The PDE2A has specificity for cAMP with Km of 2.4 μ M and its activity is modulated by the presence of cGMP. In contrast, the PDE2B does not hydrolyze cGMP nor is its activity altered by the presence of cGMP. Western blot analysis of PDE2A reveals a variety of tissues including neocortex, cerebellum, heart, kidney lungs, pulmonary artery and skeletal muscle (2). The PDE2A expression was evident in venous arterial endothelial cells but not in arterial endothelial cells (2). The PDE2A expression was also noted in corpus cavernosum along with a wide repertoire of other PDEs enzymes (3).

FabGennix provides antibodies to many PDE family members including PDE-selective, family subtype-selective and family-subtypevariant selective antibodies for the detailed analysis of cyclic nucleotide signaling pathways. We have generated two PDE2Aselective antibodies that were generated from epitopes close to N- and C-terminal regions using cyclic peptide methodology that results in higher titer and specificity (4). The polyclonal antibody (PD2A-101AP) labels a 103 kDa PDE2A protein in PDE2A positive control (PC-PD2A). The PDE2A-specific antiserum has no cross reactivity against PDE2B or with other PDE family members tested so far. Antibodies can be conjugated to secondary enzymes or fluorophores upon request at nominal costs. Antigenic blocking peptides are also available for immunocompetition/immunodepletion assays. For a complete listing of all FabGennix products and services please visit http://fabgennix.com.

References

- 1. Rascon A., Soderling SH, Schaefer and Beavo JA. Proc. Natl. Acad. Sci. USA 2002; 299, 4714-4719.
- 2. Sadhu K., Hanseley, K, V. L. Florio, and Folda, SA. J. Histochem and Cytochem. 47, 895-906, 1999.
- 3. Kuthe A., Widenroth, A., Magret J. H., et. al., J. Urol. 2001. 165, 280-283.
- 4. Farooqui S. M. Hamdi A., Brock J., Prasad C. J. Neurochem 57;1363-369, 1991.

* For users who may require large amounts of the products listed above, please inquire about bulk material discounts. This Product is for Research Use Only and is NOT intended for use in humans or clinical diagnosis.